



*keeping
you safe*

New **EM**sharps container





New **EM**sharps container

The **EM**sharps product line is designed for collecting used disposable syringes and other medical supplies. It has an easy opening for everyday usage and a permanent closure which cannot be opened after the products' final close.

- » Puncture and leak resistant
- » Shock resistant
- » Temporary closure
- » Permanent closure



The **EM**sharps are 100% PP and can be autoclaved or incinerated. They burn clean, emitting only carbon dioxide and water vapors.





» Body / Lid Stacking

The containers can be stacked in order to optimize transportation before use. When ready, the units can be easily mounted.

» Temporary Closure/ Final Closure

The **EM**sharps products offer temporary closure for user safety and to keep contents secure when not in use. Final closure ensures optimal sharps containment after use and during transportation, without the need for additional packaging.



» Usage open/ needle remover

The closure system is designed to assure the safe removal of needles and blades.



» Full opening lid



EMsharps 10 Lt



EMsharps 5 Lt



EMsharps 3 Lt



EMsharps 1 Lt

	Dim. L x W x H	Total Capacity	Useful Capacity	Max. Load	Box L x W x H	Euro Pallet 1200x800mm
EMsharps 1L ref. 7010	145 x 115 x 150 mm	1000 ml	750 ml	0.7 Kg	800x300x550mm 128 uni.	1200x800x1200mm 8 Box 1024 uni.
EMsharps 3L ref. 7030	195 x 130 x 235 mm	3000 ml	2250 ml	2 Kg	800x400x550mm 84 uni.	1200x1000x1200mm 6 Box 504 uni.
EMsharps 5L ref. 7050	280 x 155 x 215 mm	5000 ml	3750 ml	3 Kg	800x600x550mm 80 uni.	1200x1000x1200mm 4 Box 320 uni.
EMsharps 10L ref. 7100	280 x 195 x 305 mm	10000 ml	7500 ml	5 Kg	800x600x550mm 48 uni.	1200x1000x1200mm 4 Box 192 uni.

***EMsharps** is developed according to standard BS 7320:1990
All specifications subject to revision without notice.

EMsharps container

a product by



Please contact us for more
information and technical details.

emsharps@embalnor.pt

t. +351 244 800 200

f. +351 244 800 205